
## GB Patient Instructions

These exercises have been prescribed by a qualified eye care professional to treat a problem with your convergence. Regular and correct practice should improve your convergence to a normal range within six weeks, for most people. During and immediately after practicing it is common for the eyes to feel uncomfortable with the effort of doing the exercises, and a mild frontal headache may be experienced. These are positive signs

that you are exercising the eye muscles well.

Practice little and often for best results. No more than 2-3 minutes at a time, at least twice a day is recommended. Always relax the eyes after practicing by closing them or by looking into the distance for a few seconds.

**Dots:** First slightly bend the card lengthways along the central crease so it forms a wide V shape. Hold the card so the short end touches the end of the nose and the card is kept central and tilted slightly downwards. Look along the line at the farthest dot, you should see one dot and the line now appears to be a V with the dot at the bottom of the V. Now slowly look along the line to each of the other dots in turn. If your eyes are converging correctly you will see each dot at the cross over point of an X made by the line.

If the dot you are looking at is double or appears very wobbly, try hard to bring the two dots together, go back one dot further away and try again a few times. In time and with regular practice it will be possible to see more and more dots correctly. You have achieved normal convergence when you can correctly look at each dot starting at the end and finishing with the dot closest to the nose and hold each one to the count of five.

**Letters:** You can use letters on the reverse in the same way as dots. Turn the card over and hold it so the letters are in line with your nose. Both rows of letters can be used by turning the card so the letters are seen the correct way up.

The user is required to report any serious incident that has occurred in relation to this medical device to the manufacturer and to the national Competent Authority.

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