Instruction Manual

Automatic

Upper Arm Blood Pressure Monitor

Medical Disclaimer

manual and product are not meant as a substitute for advice provided by your doctor.

Model No. HI 858DK

You are not to use the information contained herein, or this product for diagnosing or treating a health problem or prescribing any medication. If you have or suspect that you have a medical problem, promptly consult your healthcare provider.

Intended Use

This device uses the oscillometric method to automatically measure systolic and diastolic blood pressure as well as heart rate

The measurement position is at human being's arm.

All values can be read out in LCD panel, and pulse waveform detected by AFib detection area will

The device is designed for home use and recommended for use by adults aged 18 years and older with upper arm circumference ranging from 9 ~ 13" (approx. 23 ~ 33 cm)

This device also equipped an AFib detection feature to collect and analyse pulse, if the specific characteristics of atrial fibrillation (AFib) is presenting during the measuring, the device will give a warning signal with the reading.

About Blood Pressure A. What is blood pressure?

A. what is blood pressure is the measurement of the force of blood pushing against the walls of the arteries. Arterial blood pressure is constantly fluctuating during the course of the cardiac cycle. The highest pressure in the cycle is called the systolic blood pressure, and represents the pressure in the artery pressure in the cycle is called its system bodo pressure, and opressing and pressure in the art of the system of t

Many factors such as physical activity, anxiety or the time of day, can influence your blood pressure. Blood pressure is typically low in the mornings and increases from the afternoon to the evening. It is on average lower in the summer and higher in the winter.

B. Why is it useful to measure blood pressure at home? Having one's blood pressure measured by a doctor in a hospital or a clinic, is often associated with a phenomenon called "White Coat Hypertension" where the patient becomes nervous or anxious, thus raising his blood pressure. There are also numerous other factors that might cause your blood The raising instruction by a provide the second state of the raise of the raise of the raise of the raised at a specific time of day. This is why medical practitioners recommend home monitoring as it is important to get readings of blood pressure during different times of the day to really get an idea of your real blood pressure.

Medical practitioners generally recommend the "Rule of 3", where you are encouraged to take your blood pressure three times in a row (at 3 ~ 5 minute interval), three times a day for three days. After three days you can average all the results and this will give you an accurate idea of what your blood pressure really is

A. WHO blood pressure classifications: Standards for assessment

of high or low blood pressure without regard to age, have been established by the World Health Organization (WHO), as shown in the chart. However this chart is not exact for classification of blood pressure and it's intered to be used as a guide in understanding non-invasive blood pressure measurements. Please consult with your physician for proper diagnosis.

B. Variations in blood pressure:

Individual blood pressures vary greatly both on a daily and a seasonal basis. These variations are even more pronounced in hyper tense patients. Normally the blood pressure rises while at work and is at its lowest

blood pressure rises write at work and is a no of during sleeping period. (hyper tense: means a person who has high blood pressure symptom.)

The graph below illustrated the variations in blood pressure over a whole day with measurement taken every five minutes.

The thick line represents sleep. The rise in blood pressure at 4 PM (A in the graph) and 12 AM (B in the graph) correspond to an attack of pain

1 Mar M

MMM

~Lm

About Atrial Fibrillation A. What is Atrial Fibrillation (AFib)^[1]?

A. What is Atrial Fibriliation (AFID)²² f Atrial fibriliation is the most common type of arrhythmia. An arrhythmia is a problem with the rate or rhythm of the heartbeat. During an arrhythmia, the heart can beat too fast, too slow, or with an irregular rhythm. A fibriliate. The term "fibriliate" detection is ginals cause the heart's two upper chambers—called the atria. —to fibriliate. The term "fibriliate" means to contract very fast and irregularly. In AF, blood pools in the atria. It isn't pumped completely into the heart's two lower chambers, called the in AF, blood pools in the atria. It isn't pumped completely into the heart's two lower chambers, called the in AF.

m - π_c, υπουρ μουs m use arra, it isn't pumped completely into the heart's two lower chambers, called the ventricles. As a result, the heart's upper and lower chambers don't work together as they should. People who have AF may not feel symptoms. However, even when AF isn't noticed, it can increase the risk of stroke. In some people, AF can cause chest pain or heart failure, especially if the heart thythm is very rapid. AF may happen rarely or every now and then, or it may become an ongoing or long-term heart problem that lasts for years.

B. How does AFib impact my family or me? One in every six strokes is AFib-related. Whilst individuals above the age of 65 are more likely to have AFib, individuals as young as 40 can exhibit AFib. Early diagnosis can help reduce the risk of a stroke. Knowing your blood pressure and knowing whether you or your family members have AFib can help reduce the risk of stroke. HL858DK AFib detection provides a convenient way to detect the specific characteristics for AFib during the measuring pulse period

C. Risk factors you can control: High blood pressure and AFib are both considered «controllable» risk factors for strokes. Knowing your blood pressure and knowing whether you have AFib is the first step in proactive stroke prevention.

[1] : National Institutes of Health, U.S. Department of Health and Human Services.

Precautions

* Do not use this manual and product as a substitute for advice, diagnosing or treating a health problem or prescribing any medication by your doctor. If you have a medical problem, promptly consult your healthcare provider. Read the Instruction Market and the problem, promptly Read the Instruction Manual thoroughly before measuring and keep it at hand for your reference.

at any time. at any time. This device does the oscillation of the source of the sourc

necessarily by a physician or operator.
 This monitor is not intended for use in the MR environment.

 \odot Do not take a measurement in a low (less than 41 °F/5 °C) and high (more than 104 °F/40 °C) temperature, nor in a place outside humidity ranges (15 % \sim 93 % R.H.), and atmospheric pressure ranges (700 \sim 1060 hPa), or you may get inaccurate readings. \odot Wait 30 \sim 45 minutes before measurement if you've just consumed caffeinated beverages or smoked cigarettes

 \odot Rest at least 5 ~ 10 minutes before taking a measurement. \odot To allow your blood vessels to return to the condition prior to taking the measurement, please wait at least 3 ~ 5 minutes in between measurements. You may need to adjust the wait time

according to your personal physiological situation. • We recommend you using the same arm (preferably the left arm) and measuring around the

Sit down comfortably and place your elbow on the table with your feet flat on the floor. Please

 of a dwin consist and a dwing measurements.
 O Keep the cuff at heart level. Relax your hand with the palm facing up.
 O Perform measurements in a quiet and relaxed environment at room temperature.
 O Do not move or shake the device during a measurement. Please keep quiet and do not talk during measurements.

• This product is not suitable for:

Pregnant women The diagnosis of arrhythmia Undergoing intravenous injection on any limb Currently in a dialysis treatment In pre-eclampsia condition

• For those who have had a mastectomy or lymph node clearance, it is recommended to take a When used among medical electronic equipments on the same limb, pressurization of the cuff

may cause temporarily malfunction to other devices.

 ⊙ Keep in mind that blood pressure naturally varies from time to time throughout the day and is affected by lots of different factors such as stress, eating, smoking, alcohol consumption, medication, and physical activity, etc.

mencication, and physical activity, etc. © Normally the blood pressure rises while at work and is at its lowest during sleeping period. © Blood pressure measurements should be interpreted by a physician or a trained health professional who is familiar with your medical history. Using the unit and recording the results

regularly for your physician to interpret, you will keep your physician informed of the continuing changes in your blood pressure \odot If you have one of the circulatory problems as arteriosclerosis, diabetes, liver disease, kidney

disease, severe hypertension, peripheral circulation....., please consult your healthcare professional before using the device Results are not intended for direct diagnosis. Please consult with a physician if you have any

Questions are not interface of direct diagnosis. Tease consult with a physician in you have any questions or concerns about your results. © Blood pressure measurements determined with this device are quivalent to those obtained by a


trained observer using the cuff/stethoscope auscultation method and are within the accuracy limits prescribed by the American National Standard for Manual, electronic, or Automated Sphygmomanometers. If the cuff is worn incorrectly, or the shape of the upper arm is special (for example, the

In the curit is work incorreculy, of the starge of the upper aim is special for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive gap might occur between the arm cuff and the arm, and it might lead to measurement errors or inaccuracies. If you have any question about the condition of cuff wearing and/or measurement result, please consult your healthcare professional. • The applied part is cuff.


*Attention

 Do not use the device on infants, children, or those who cannot express their own intention. To avoid accidental strangulation, keep this product away from children and do not drape tube around neck. 2. The medical device should not used adjacent to or starked with other equipment. In case adjacent or stacked with other equipment. In case adjacent or stacked with other equipment. In case adjacent or stacked use is necessary. The medical device should be observed to verify normal operation in the configuration in which it will be used.

configuration in which it will be used. 3. Consider the electromagnetic compatibility of the device (ex. power disturbance, radio frequency interference etc.) Please use it indoor only. 4. Over high frequency measurements may result in blood flow interference, which is likely to cause uncomfortable sensations, such as partial subcutaneous hemorrhage, or temporary numbness to your arm. In general, these symptoms should not last long. However, if you do not recover in time, please seek your medical practitioners for help.



*Caution! Substitution of a component different from that supplied might result in measurement error



	Low Ballery	polarity reverses.
Symbol		→ We suggest you replace all batteries with new ones, and make sure the +/- polarities are properly positioned.
		Once pulse is detected, the symbol flashes with each pulse beat.
	Pulse Symbol	→ Our suggestion: Please do not talk or move during measurements.
AVG.3	Memory Average	This symbol appears when LCD displays average value of last 3 readings.
Irregular Heartbeat Detector	Irregular	This symbol appears for 1.5 minutes when the user was talking, moving, shaking, or an irregular heart beat was detected during measurements.
		→ Our suggestion: Please do not talk or move during measurements. Repeat the measurement after resting for at least 5 minutes, and restart your measurement while sitting down comfortably and quietly.
	BP Category Indicator	The arrowhead points out the specific BP Category that your measurement reading fits in.
	Excessive Body Motion Detector	Displayed if body movement is detected during measurement, especially, the movement on the arm the blod pressure monitor is worn on. Besides, if colf is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of the forearm), excessive age might exist between the arm culf and the arm. Notice: The measured blod pressure reading may not be accurate if the icon is displayed.
	Cuff Fitness detection Symbol	Displayed if the cuff was wrapped incorrectly, which is too tight or too loose. This is the function aid in detecting if the cuff is wrapped properly.

This symbol appears when there is no detection of Atrial Fibrillation $\mathcal{A}\mathcal{A}$ eature elated Symbol This symbol appears when there is detection of Atrial Fibrillation

Features

BP Category Indicator

This device is equipped with BP Category Indicator which classifies your blood pressure measurements into six stages (Optimal to Severe hypertension) as shown in below chart:

Sta Pre	ages of Blood essure Levels	Systolic (mmHg)	Diastolic (mmHg)	Color	Recommendations by SIGN 49: Hypertension in older people
Grade 3	Severe Hypertension	≥180	≥110	Red	Confirm immediately and repeat BP in one day and again within one week depending on clinical situation.
Grade 2	Moderate Hypertension	160~179	100~109	Red	Serial blood pressures repeated within one month.
Grade 1	Mild Hypertension	140~159	90~99	Red	Provide advice about lifestyle modification and confirm within two months.
Н	igh-Normal	130~139	85~89	Orange	Provide advice about lifestyle modification and recheck in one year.
	Normal	120~129	80~84	Yellow	Recheck in 2 - 5 years.
	Optimal	<120	<80	Green	(patients aged > 75 years offered annual health check)

After each measurement is completed, LCD display will show your position automatically on the six segments of the bar indicator which corresponds to BP Category Indicator.

IASTOLIC (mmHg) 110	Seve	ere hyp	perter	sion				- [- • E- Red
105	Moderat	te hyp	ertens	ion					— ▶ ■ — Red — ▶ ■ — Red
100 95	Mild hypert	oncior							- Ned
90		ension							─ ► Orange ─ ► Yellow
85	High normal								• I =
80	Normal — Optimal —				_				- Green
	120	130	140	150	160	170	180 S	YSTOLIC (mmHa)	

*Note!

When a person's systolic and diastolic pressures fall into different categories, the higher category should apply.

appy. e.g. systolic pressure 181 & diastolic pressure 99 → Red category (Severe Hypertension) e.g. systolic pressure 110 & diastolic pressure 95 → Red category (Mild Hypertension)

*Note!

The above table is not exact for classification of blood pressure and it's intended to be used as a guide in The above table is not exact for classification or blood pressure and its intended to be used as a guide in understanding non-invasive blood pressure measurements. Usually this is not a cause for concern; however we recommend you consult with your physician for proper diagnosis or seek medical advice according to our recommendation mentioned above. Please note that the device does not appropriate to diagnose hypertension, and it is only for user reference on blood pressure

Irregular Heartbeat Detector

The symbol will appear on screen indicating a certain heartbeat irregularity was detected during measurement.

The heartbeat rhythm that is more than or less than 25% from the average rhythm is usually defined as an irregular heartheat rhythm

Talking, moving, shaking or an irregular pulse during the measurement can result in the appearance of this symbol

Usually this is not a cause for concern, however if the symbol appears often, we recommend vou seek medical advice.

And please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

*Note!

• The pulse display is not suitable for checking the frequency of heart pacemarkers. If a certain pulse Ine pulse display is not suitable for checking the trequency of heart pacemarkers, if a certain pulse
 irregularity is detected during measurement often, we recommend you seek medicail advice
 As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature bats
 and atrial fibrillation or any other special conditions you should check with your physician before using your

device. • The IHS function is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmic problem. In order to filter the unstable status of user and avoid affecting the detection of heart and promote provide the method of a state of the second of

• At least 3 beats with at least 25% difference from the average heart beat interval will generate the IHB icon

Atrial Fibrillation Detection Feature

This device is equipped an AFib detection feature to collect and analyse pulse signal frequency from the user's finger, it provides a convenient way to screen AFib condition during measurement. The detection of atrial fibrillation is determined by the collection from a period of pulse signals from the AFib detection area

If atrial fibrillation is detected during measurement, the AFib symbol displayed. If AFib symbol appears more frequently, we recommend the patient to seek professional medical advice

This device does not replace a cardiac examination, but attempt to detect atrial fibrillation at the early stage.

*Note!

 Sometimes the device will detect atrial fibrillation even when it is not there. This could Owneamest and evence while detect at the instance even where it is to there it must be that and and finger move during the measuring or another rhythm problem is present. Keep the hand and finger still during the measuring.
This device may not detect attrait fibrillation in people with pacemakers or

defibrillators. • We recommend you consult with your physician for proper diagnosis or seek medical advice according to our recommendation mentioned above. Please note that the device does not appropriate to diagnose atrial fibrillation, and it is only for user reference on blood pressure monitoring.

Installing Batteries

When LOW BATTERY SYMBOL D appears on the display, or no reaction toward operation, please change batteries.

Replace all worn-out batteries with new ones and do not mix new and used batteries. Do not mix alkaline, standard (carbon-zinc) or rechargeable (cadmium) batteries either. Such action may shorten the battery life or cause the device to malfunction

Slide the battery cover and insert 4 AA (I R6) alkaline batteries into the battery compartment as shown on the figure below. Make sure the polarities "+" and "--" ends are coinciding with similar markings engraved on the battery housing.

*Attention

*Note!

*Note!

skin

vour heart

*Note!

cleared)

12).

standby mode

User 1 or 2

on the screen

process is completed.

VOU

120

9

D

HIM

damage of the unit.

Rating: Input: 100-240V, 50/60 Hz, 0.2 A Output: 5V, DC, 1A, O

to disconnect the adapter plug.

the inside of the cuff.

above the inner elbow.

Applying the Cuff

Using the AC adapter

No batteries are needed when operating with an AC adapter.

Recommend Adapter specification, do not use otherwise: Model: FranMar International, FRM06-S05-EU

- Attention!
 Batteries are hazardous waste. Do not dispose of them together with the household garbage. Please discard worn-out batteries to the recycling site according to local regulations.
 Keep the battery away from children in case they choke on it.
 To prolong the battery life and prevent damage caused by leakage, remove the batteries from the device if the device is not to be used for a long period.
 The device will keep the last measuring results after changing batteries, please reset date and time.
 Please replace all worn-out batteries with new ones when you are operating the Atrial Fibrillation Detection feature, and the LOW BATTERY SYMBOL → appears on the display.

This monitor is designed for operation with batteries or an AC adapter. Please use only a compatible AC adapter with required voltage and current as indicated

Please unload the batteries when operating with an AC adapter for an extended period of time.
 Leaving the batteries in the compartment for a long time may cause leakage, which may lead to

When you use the blood pressure monitor with AC adapter, do no position the device to make it difficul

• Wrap the cuff on a bare arm or over thin clothing. Thick clothing or a rolled up sleeve

incorrect measurement result. \odot Press your brachial artery approximately 1 inch (2 ~ 3 cm) above the elbow on the

Fit the cuff snugly, leaving enough space for 1 inch (2 ~ 3 cm) between the inner elbow and the lower edge of the cuff, or the measurement may not be accurate.
 This monitor comes with one size of arm cuff. 9" ~ 13" (23 ~ 33 cm).

This monitor comes with one size of arm cuff. 9" ~ 13" (23 ~ 33 cm).
In case the cuff kept pumping up non-stop unwrap the cuff at once.
Do not wrap the cuff around any body part other than your arm.
The device is not supposed to be used when your arm is wounded or injured.
If you have any infectious skin disease or the device is used by users with infectious skin disease, please do not continue using the device.
Before using the device, user should check the appearance of cuff. If you notice blood or other soil on cuff, please do not use this device.
If there is no of above situations, please dispose the device without reuse.
Do not use this device if your wrist (Arm) has any wound or injury, especially after surgery on the wrist (Arm). Otherwise, it may cause infection at the surgical site. Please use the device after the wound has healed.

End of Cuff

Arrow

D-rina

00

Press Start

nad

Press Star

B

15.00

STANDBY MOD

no.

Press Start

*888 *388

m

Slide the end of arm to determine where your strongest pulse is.
 Slide the end of arm cuff furthest from the tube through

lie over the brachial artery on the inner part of the arm. The bottom edge of the cuff should be 2 ~ 3 cm (approx. 1 inch)

 Pull the end of the cuff so that it tightens evenly around your arm, allow room for 2 fingers to fit between the cuff and your arm.
 Please make sure the cuff do not slip during measurement, and the arrow falls within the Proper Fit Range.

Sit on a chair, back and arm supported, and lay your forearm on the table so that the cuff is at the same level as

outside of the cuff and metal ring will not touch your Proper Fit Range

the metal ring to a loop. The smooth cloth should be on

If the cuff is located correctly, the velcro will be on the

⊙ Put left arm through the cuff loop. The tube should

When the cuff is positioned properly, press the velcro firmly against the pile side of the cuff.

Relax your arm and turn your arm upward.

• Make sure there are no kinks in the air tube

Measurement Procedure

A Put in 4 AA "LR6" (1.5 V) alkaline batteries

Setting Year, Date and Time

Taking a Measurement

B. All segments appear on the screen for 3 seconds.

and MINUTE (00,01.....,59) by following Step B.

A. Before measurement, press USER button to select

B-1. Start a Measurement (With the blood pressure measurement and AFib detection feature):

1. With the cuff wrapped around your upper arm, and place the finger of the opposite hand (index finger

recommended) gently on the AFib detection area, then

START button to start measurement. All display units appear

2. After all symbols disappear, the display will show "00". The monitor is

3 As the cuff inflates, the monitor automatically determines your ideal

inflation level. This monitor detects your blood pressure and pulse rate during inflation. The Heartbeat Symbol (•) flashes at every heartbeat

detection area. Remain still and do not move until the entire measurement

monitor also detects your pulse signals by the Advanced IHB

"Ready to Measure" and will automatically inflate to the level that is right for

make sure the finger place at the correct position, press

C. The monitor will automatically turn to sleeping mode (all LCD segment

A.Press () button ("YEAR" flashes). Press + USER button to adjust YEAR value

B. Press () button ("MONTH" flashes). Use + USER button to adjust MONTH (1, 2, 3,.....

C.Continue to set current DATE (varies from 1 to 31), HOUR (1, 2.....12PM, 1PM.....,12)

D.When settings are done, press (button to confirm the settings. The device turns to

Switch on the Monitor

LCD screen displays your systolic rate, diastolic rate, pulse, BP Category Indicator, AFib Detection feature-related symbol and Irregular Heartbeat Detector symbol (if any) 4 with date and time

*Notal

- Do not inflate the cuff until it is wrapped around your upper arm.
 Please clean the finger and make sure nothing is covered on it before taking the Atrial Fibrillation Detection measurement.
- The Atrial Fibrillation Detection measurement is not supposed to be used when your finger is wounded If the cuff does not stop inflating, remove the cuff at once.
- Please do not move and wait for the blood pressure measurement, as well as the AFib detection feature results displayed.
- 5-1. If AFib is not detected during the measurement as normal result, the 120 LCD display will be as below

へべん 5-2. If AFib is detected during the measurement, the "AFib" symbol will be displayed on the LCD as below: 12:01

6. The blood pressure measurment is completed and without any operation for 1.5 minutes, device automatically shuts off

B-2. Start a Measurement (Only with the AFib detection feature):

- When press the START STOP key, the measurement of blood pressure and AFib detection will be both activated. If user just want to take an AFib detection measurement, please place the finger gently on the AFib detection area, then make sure the finger cover the AFib detection area, press the START key again.
- 2. The display will show "User 1 or 2 with date and time". The monitor automatically detects your pulse signal. Remain still and do not no i move your finger until the entire measurement process is completed
- When the measurement is completed, LCD screen displays AFib 3 Detection feature-related symbol with date and time. $\overline{}$

*Note! If user only takes the Atrial Fibrillation Detection measurement, data can Not be stored.

Memory Function

Storing data

The monitor can store up to 240 memories for 2 users, and automatically replace the oldest data with new one

The	Condition of	AFib detection feature					
Mea	surement	Success	Error				
The blood pressure measurement	Success	This condition will be stored in the memory and recalling data will display as below:	E4 Try Again E4 Try Again This condition will be stored in the memory, however recalling data will display as below:				
The	Error	This condition will NOT be stored in the memory.					

Recalling data A. Press USER button to select User 1 or 2. B. Press M button to enter Memory Mode. LCD displays average of last 3 measuring results first. C. Press M button again, LCD displays the latest measuring result. Use M and + button to scroll through all stored measuring results. D. To stop reading memories, press STAP Mode.	No Data
Erasing data A. Press USER button to select User 1 or 2. B. Press M button to enter Memory Mode. C. Press and hold + and ④ buttons at the same time, the data will be	
erased automatically. D. To confirm the data in the selected user has been erased, press M button and no data should appear. *Note! Once deleted, your data can NOT be restored.	[1]

Storage and Maintenance

General Use Do not in any way twist the cuff.

Do not drop the product and avoid any strong impacts.

Maintenance

To ensure that your device is in optimal use and to avoid damage, please refer to the following instructions: Clean the device and cuff with a soft dry cloth, or

- Use a dry cloth with water to clean the device (not directly flush, do not soak in water, and
- Dee a di y cluti war water, and bella the device (not allecty instit, do not soak in water, and hold the device dry), or
 Do not use detergent or any strong chemicals to clean the device.
 Make sure the cuff is completely dry before using.
 According to the use environment of the sphygmomanometer, the recommended disinfection method and frequency are as follows:
- Only use it yourself (home use), it can be cleaned at ordinary times, and wipe it once a
 month with a commercially available 75% alcohol cotton sheet (for the cuff) for more than 30 seconds each time.
- 30 seconds each time. If it is used for more than one person (home use), it can be cleaned at ordinary times. It is disinfected once a week (for the cuff belt) with a commercially available 75% alcohol cotton sheet, for more than 30 seconds each time. After cleaning / disinfection/ before use, please make sure that there are no blood stains or soil on the LCD, the device and cuff, If there is any blood stains or soil, please dispose the device with each time.
- the device without reuse.
 - If it is used in a complex environment (such as a hospital) or after multiple people (non-family), please discard the old cuff and replace it with a new one.

Storage

BŪ

120

80

 $\overline{\mathcal{M}}$

80

80

- If the device is not to be used for a long time, please remove the batteries from the device (leaking of battery acid can cause the device to malfunction). Always store the unit in the storage case after use. It is intended to be transported or
- stored in a carrying case between uses. Do not place the device directly under sunlight, in high temperature, or in humid or dusty
- places

Troubleshooting

SYMBOLS/ SYMPTOMS	CONDITIONS/CAUSES	INDICATION/CORRECTION	
Unit does not turn on when START STOP button is pushed.	Worn-out batteries.	Replace them with 4 new AA (LR6) alkaline batteries.	
STOP button is pushed.	Battery polarities have been positioned incorrectly.	Re-insert the batteries in the correct positions.	
66	Cuff has been placed incorrectly.	Wrap the cuff properly so that it is positioned correctly.	
Measuring Error Symbol	Did you talk or move during measurement?	Measure again. Keep arm steady during	
appears when blood pressure value displayed s excessively low or high.	Shaking of the arm with the cuff on.	measurement.	
Measuring Error Symbol	Air circuit abnormality. Cuff tube may not be plugged into monitor correctly.	Check cuff connection. Measure again.	
Measuring Error Symbol	Inflation pressure exceeding 300 mmHg.	Switch the unit off, then measure again.	
Measuring Error Symbol	Can't determine blood pressure measurement data.	Wrap the cuff properly and keep steady. Measure again.	
Excessive Body Motion Detector	Body movement during measurement, especially, the movement on the arm the blood pressure monitor is worn on. e.g. Talking, moving or shaking of the arm with the cuff on while measurement.	Measure again. Keep arm steady during measurement.	
Notice: The measured blood pressure reading may not be accurate if	Cuff is worn improperly, or the shape of the upper arm is unusual (for example, the circumference of the upper arm differs largely from the circumference of	Wrap the cuff properly and keep steady. Measure again. If you have any question about the cuff	
the icon is displayed.	the forearm), excessive gap might be exist between the arm cuff and the arm.	wearing and/or measurement result, please consult your healthcare professional.	
Cuff Fitness detection Symbol	The cuff was wrapped incorrectly (for example too loosely or too tightly).	Please reference "applying the Cuff " section to wrap the cuff correctly	
	Finger hasn't be placed on AFib detection area.	Keep finger gently place and well-covered on AFib detection area and measure again.	
E4 Try Again Measuring Error Symbo	Finger moved away from the detection area when the measurement has not been completed yet.	Measure again and don't move away your finger before measurement completed.	
E5 Try Again	Finger press too hard on AFib detection area.	Measure again. Gently place finger on AFib detection area.	
Measuring Error Symbo	Cold Finger and weak pulse signals can't determine AFib measurement data.	Measure again. Keep finger warm and gently place on AFib detection area.	
E6 Try Again Measuring Error Symbo	Pulse signals could not be detected continuously by the AFib detection area for a period.	Place the finger on AFib detection area and keep steady. Measure again.	
	 display, just return the device to your local distributo 	I	

Limited Warranty

Warranty For Two Years from the manufacturing date

Please note that this warranty does not cover damage caused by misuse or abuse; accident: the attachment of any unauthorized accessory: alteration to the product: improper installation; unauthorized repairs or modifications; improper use of electrical/ power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair, or any other conditions whatsoever that are beyond the control of importers or distributors.

In case it is needed to have the device checked for calibration, please consult the distributor. This is recommended to be considered every two years.

Specifications

Model Number	HL858DK		
Measurement Method	Oscillometric		
Measurement Range	Pressure: 0 ~ 300 mmHg		
	Pulse: 40 ~ 199 Beats/Minute		
Accuracy	Pressure: ± 3 mmHg		
	Pulse: ± 5 % Max.		
Inflation	Automatic Inflation (Air Pump)		
Deflation	Automatic Air Release Control Valve		
Display	Liquid Crystal Display		
Memory	240 Memory Total for 2 Users		
Unit Dimensions	6.08 x 4.24 x 2.37 inch (L x W x H)		
	154.5 x 107.8 x 60.4 mm (L x W x H)		
Unit Weight	295.8 g ± 5 g (10.43 oz ± 0.17 oz)		
	(Cuff and batteries excluded)		
Cuff Size	NC-01: Normal size cuff 9 ~ 13 inch (23 ~ 33 cm)		
	UC-01: Universal size cuff 9~17 inch (23 ~ 43 cm)		
	Optional		
Storage/	Temperature: -25 °C ~ 70 °C (-13 °F ~ 158 °F)		
Transportation Environment	Humidity: ≤ 93 % R.H.		
Operation Environment	Temperature: 5 °C ~ 40 °C (41 °F ~ 104 °F)		
	Humidity: 15 % ~ 93 % R.H.		
	Atmospheric pressure: 700hPa ~ 1060hPa		
Power Supply	AA "LR6" (1.5 V) alkaline battery x 4		
	5 V 1A AC Adapter (Optional)		
Battery Life	Approx. 200 Measurements		
Product Life	5 Years (4 times per day)		
Sleeping Mode	Without any operation for 1.5 minutes,		
	device automatically shuts off.		
Accessories	4 AA (LR6) Alkaline Batteries, Arm Cuff with Tube,		
	Instruction Manual, Storage Pouch		

Note

artially applied) NI SO 81060-2 Non-invasive sphygmomanometers - Part 2: Clinical investigation of intermittent automated measurement type CC 80601-2-30 Medical electrical equipment - Part 2-30: Particular requirements for basic safety and essential performance of utomated non-invasive sphygmomanometers.
phygmomanometers. N ISO 81060-1 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type
N 1060-4 Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive
Verformance standards: N 1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring vertems.
N 60601-1-2 Medical electrical equipment part 1-2: General requirements for safety- Collateral standard: Electromagnetic ompatibility- Requirements and tests
Safety standard: N 60601-1 Medical electrical equipment part 1: General requirements for safety And essential performance. MC standard:
This blood pressure monitor complies with the EC Directive (93/42/EEC) and bears the ZE mark. This blood pressure monitor also complies with mainly following standards CE 0197

Explanation of symbols :

Symbol	Explanation	Health & Life Information
ĆE	CE conformity marking	-
0197	Notified Body (NB) number	-
8	Refer to instruction manual/ booklet	-
Ŕ	TYPE BF Applied Part	-
Ć	To avoid inaccurate results caused by electromagnetic interference	Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the device. Otherwise, degradationof the performance of this equipment could result.
X	Waste of electrical and electronic equipment (WEEE)	Discard the used product to the recycling collection point according to local regulations
	Manufacturer	HEALTH & LIFE Co., Ltd. 9F, No.186, Jian Yi Road, Zhonghe District 23553, New Taipei City Taiwan www.healthandlife.com.tw
\sim	Date of Manufacture	м үүүү-мм
EC REP	Authorized representative in the European Community	EC REP EMERGO EUROPE Prinsessegracht 20, 2514 AP The Hague, The Netherlands
SN	Serial number	SN YYMMXXXXXX
IP22	Ingress Protection Rating	First characteristic numeral-Degree of protection against access hazardous parts and against sold foreign objects N1=2 (Protecte against solid foreign objects of 12.5 mm Ø and greater) Second characteristic numeral-Degree of protection against ingress of water N2=2 (Protected against vertically failing water drops when ENCLOSURE tilled up to 15')
^{93 %}	Humidity limitation (Storage/Transportation condition)	R.H.: ≦93 %
-25 °C (-13 °F)	Temperature limit (Storage/Transportation condition)	Temperature: -25 °C ~ 70 °C (-13 °F ~ 158 °F
700 HPa	Atmospheric pressure limitation (Operating condition)	Atmospheric pressure: 700 hPa~1060 hPa
^{93 ۹}	Humidity limitation(Operating condition)	R.H.: 15 % ~ 93 %
5 °C (41 °F)	Temperature limit (Operating condition)	Temperature: 5 °C ~ 40 °C(41 °F ~ 104 °F)
$((\cdot,\cdot))$	Non-ionizing electromagnetic radiation	-

Device information:

- Internally powered equipment - Not suitable for use in presence of flammable anesthetic mixture with air or with Oxygen or nitrous oxide - Continuous operation with short-time loading

HEALTH & LIFE CO., LTD.

9F, No. 186, Jian Yi Road, Zhonghe District 23553, New Taipei City, Taiwan www.healthandlife.com.tw

Appendix

		the electron	anotio or in	opponto ¹¹ -1	ad below, and abould only be used in such any			
Emissions test		Complian			ed below, and should only be used in such environments:			
		Compliance			omagnetic environment-guidance			
RF emissions CISPR 11		Group 1		Therefore,	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.			
RF emissions CISPR 11		Class B		The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply				
Harmonic emissions IEC 61000-3-2 Voltage fluctuations/		Class A		network that supplies buildings used for domestic purposes.				
flicker emissions IEC 61000-3-3		Complies						
Guidance and	manu	facturer's	s declar	ration –	electromagnetic immunity			
The device is intended	l for use in	the electroma	ignetic envir	onments lis	ed below, and should only be used in such environments:			
Immunity test		60601 t level	Compl lev		Electromagnetic environment-guidance			
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV co discharge ± 15 kV a discharge) ir	± 8 kV con discharge ± 15 kV air discharge		In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15 °C-35 °C, Relative Humidity: 30 %~60%.			
Power frequency (50 or 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz		30 A/m 50 or 60 Hz		Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.			
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output		± 2 kV for power supply lines ± 1 kV for input/output		Mains power quality should be that of a typical commercial or hospital environment.			
Surge IEC 61000-4-5	AC Power port		AC Power port ±1 KV Line to Line		Mains power quality should be that of a typical commercial or hospital environment.			
interruptions and voltage variations on power supply input lines	0% UT; 0.5 cycle At 0°,45°,90°,135°, 180°,225°,270°and 315°.		0% UT; 0.5 cycle At 0°,45°,90°,135°, 180°,225°,270°and 315°.		Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruphons, it is recommended that the device be powered from an uninterruptible power supply or a battery.			
IEC 61000-4-11		% UT; 1 cycles 0 % UT; 25/30 cycles		cycles 5 cycles	uninterruptible power supply or a battery.			
	0%117.2	T; 250/300 cycle 0 % UT; 2		0 cvcle				
Conducted RF IEC 61000-4-6	3V rms At 0.15-8	3V rms			Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, thar the recommended separation distance calculated from the			
	6V rms At ISM & Amateur	Freq.	6V rms At ISM & R Amateur Fr		equation applicable to the frequency of the transmitter.			
Radiated RF IEC 61000-4-3 (Proximity fields from RF wireless communications equipment IEC 61000-4-3)	9-28V/m a 385-6000 Mode and Modulati system sh as specifi IEC60601 for proxim from RF v communic equipmer	MHz at MHz, d other on. The all be tested ed in 1-1-2 Table 9 ivtreless cations at using the ods specified	9-28V/m at 385-6000M er Mode and c he ke tested system sha as specified for proximit for mroximit for mroximit for mroximit ssc for generic equipment specified test methoc		Recommended separation distance Considering to reduce the minimum separation distance, based on RISK MANACEMENT, and using higher MMUNITY TEST LEVELS that are appropriate for the acculated using the following equation. $E = 640, F^{T}$ where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVELS in V/m. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.			
					Interference may occur in the vicinity $(((\cdot)))$ of equipment marked with the following symbol:			

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. surcures, objects and people. a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measure strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as

observed to verify normal operation. If admirinal performance is observed, additional measurements of the device. Foreining or relocating the device. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Test frequency (MHz)	Modulation	IMMUNITY TEST LEVEL (V/m)	
385	Pulse modulation 18 Hz a)	27	
450	FM ± 5 kHz deviation 1kHz sine ^b)	28	
710			
745	Pulse modulation 217 Hz a)	9	
780			
810			
870	Pulse modulation 18 Hz ^a)	28	
930			
1720			
1845	Pulse modulation 217 Hz ^a)	28	
1970			
2450	Pulse modulation 217 Hz ^a)	28	
5240			
5500	Pulse modulation 217 Hz ^a)	9	
5785			

a). The carrier shall be modulated using a 50 % duty cycle square wave signal.
b). AS an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.