EN

HemoCue® Hb 801 Microcuvettes, 4x50 Microcuvettes

HemoCue Hb 801 System consists of HemoCue Hb 801 Analyzer and HemoCue Hb 801 Microcuvettes. Read and follow this package insert together with the operating manual (provided with the analyzer), for optimum performance and safety. For more information, contact the local distributor or HemoCue AB. Only use HemoCue Hb 801 Analyzer together with HemoCue Hb 801 Microcuvettes.

Intended Use

The HemoCue Hb 801 System is intended for the quantitative determination of hemoglobin in capillary or venous whole blood in point-of-care and blood donation settings. The HemoCue Hb 801 System is for professional *in vitro* diagnostic use only.

IVD Medical Device Directive

The HemoCue Hb 801 Microcuvettes comply with the IVD Medical Device Directive 98/79/EC and carry the CE mark.

Composition

The microcuvettes are made of polystyrene plastic, and contain no active ingredients.

Precautions

Always wear protective gloves. Handle blood with care, as it may be infectious. Follow local safety procedures for disposal of used microcuvettes.

HemoCue Hb 801 Microcuvettes are for single use only.

For system information regarding precautions, read the operating manual.

Storage and Handling

Store the microcuvettes in 10–40 °C (50–104 °F). Microcuvettes in the vial (opened or unopened) are stable until expiration date, printed on the package. Use the microcuvettes prior to expiration date.

Microcuvettes in the vial (opened or unopened) can be stored for a shorter period of time (6 weeks) between -18–50 °C (0–122 °F).

Keep all unused microcuvettes in the original package.

Measuring

Operating temperature for HemoCue Hb 801 System is 10–40 °C (50–104 °F).

For proper use of the system, read the operating manual regarding measuring capillary and venous samples, measuring control material, cleaning and disinfection, principle of method, and technical specifications.

Make sure the microcuvette is completely filled before measurement.

Immediately after filling, insert the microcuvette into the analyzer to perform a measurement. Do not let more than 40 seconds pass.

Quality Control

The HemoCue Hb 801 Analyzer has an internal quality control, a self-test. For more information, read the operating manual. If an external quality control is required by local or other regulations, only use controls recommended by HemoCue.

Limitations

- Start measurement immediately after filling the microcuvette. Do not let more than 40 seconds pass.
- Following substances have not been found to interfere: acetaminophen (1324 µmol/L), ascorbic acid (342 µmol/L), conjugated bilirubin (15 mg/dL), unconjugated bilirubin (10 mg/dL), creatinine (442 µmol/L), ibuprofen (2425 µmol/L), leucocytes ($30 \times 10^9/L$), Intralipid® (125 mg/dL), salicylic acid (4.34 mmol/L), tetracycline (34 µmol/L), urea (42.9 mmol/L), uric acid (1.4 mmol/L), HbCO (10 %), MethHb (5 %), hemolysis (2 g/L), HbO₂ (≤ 50 %), warfarin (32.5 µmol/L), simvastatin (49 µmol/L), and platelets (2000×10^9 cells/L). The highest concentration or percentage tested is referred to in brackets. Interference studies have been performed according to CLSI EP07-A2.
- Normal blood pH and above, up to 8, do not interfere with the system.
- Venous samples (K₂EDTA or Li-Heparin) can be stored and transported for up to 24 hours in room temperature or refrigerator, before measuring.

Expected Values

Due to a wide range of conditions (dietary, geographical, etc.), which affect normal values, it is recommended that each laboratory establishes its own normal range. For general guidance to normal reference values applicable to most healthy adults and children, see for example Dacie and Lewis, Practical Haematology or Soldin, Pediatric Reference Intervals.

Specific Performance Characteristics

Repeatability and Within-laboratory Precision

Results summarized below was determined according to CLSI EP05-A3.

Control Level	N	$\bar{x}$ g/dL	Repeatability		Within-laboratory Precision	
			SD g/dL	CV %	SD g/dL	CV %
1	240	6.34	0.05	0.7	0.04	0.7
2	240	11.50	0.05	0.4	0.05	0.5
3	240	15.36	0.15	1.0	0.16	1.0

Accuracy

Results summarized below was determined according to CLSI EP09-A3.

Study	N	Min g/dL	Max g/dL	Regression Line	Correlation Coefficient (r)
1	102	1.3	23.6	$Y = 0.99X + 0.09$	1.00
2	100	11.8	18.4	$Y = 0.96X + 0.92$	0.92
3	239	1.5	25.4	$Y = 1.00X - 0.13$	1.00

1 = venous blood, compared to the HiCN reference method according to the International Council for Standardization in Haematology (ICSH)

2 = capillary blood, compared to the HiCN reference method according to the International Council for Standardization in Haematology (ICSH)

3 = multicenter study, venous blood, compared to HemoCue Hb 301

References

1. HemoCue Hb 801 Operating Manual
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 and references herein
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AAC Press; 7th edition, 2011

Symbols Used

	Caution		Batch code
	In Vitro Diagnostic medical device		Catalog number
	Do not reuse		Temperature limitation
	CE mark		Consult instructions for use
	Manufacturer		This product is covered by one or more patents
	Use by (year month day)		Date of manufacture

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Références

1. Manuel d'utilisation HemoCue Hb 801
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 et références y mentionnées
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACCPress; 7th edition, 2011

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- Manual de funcionamiento de HemoCue Hb 801
- Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 y las referencias que aquí se hacen
- CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
- CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
- CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
- CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
- Soldin, S. J., Pediatric Reference Intervals, AACC Press; 7th edition, 2011

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Bibliografia

- Manual de Funcionamento do HemoCue Hb 801
- Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 e respectivas referências
- CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
- CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
- CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
- CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
- Soldin, S. J., Pediatric Reference Intervals, AACCPress; 7th edition, 2011

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Referenser

1. HemoCue Hb 801 Bruksanvisning
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 och referenser häri
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACCPress; 7th edition, 2011

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Referencer

1. Brugsanvisningen til HemoCue Hb 801
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016, samt referencer heri
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACCC Press; 7th edition, 2011

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1. Bruksanvisning for HemoCue Hb 801
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 og referanser i denne
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACC Press; 7th edition, 2011

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3 = studio multicentrico, sangue venoso, comparato a HemoCue Hb 301

Bibliografia

1. Manuale d'uso HemoCue Hb 801
2. Dacie e Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 e relativa bibliografia
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACCPress; 7th edition, 2011

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Literatur

1. HemoCue Hb 801 Bedienungsanleitung
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12. Auflage, 2016, sowie darin enthaltene Quellenverweise
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACC Press; 7. Auflage, 2011

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Referenties

1. Gebruikershandleiding HemoCue Hb 801
2. Dacie and Lewis, Practical Haematology, Elsevier Limited, 12th Edition, 2016 en referenties hierin
3. CLSI H15-A3 Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition
4. CLSI EP07-A2; Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition
5. CLSI EP05-A3; Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition
6. CLSI EP09-A3; Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition
7. Soldin, S. J., Pediatric Reference Intervals, AACC Press; 7th edition, 2011

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