

# DC 430 M A S

## Dual Frequency Body Composition Monitor with Integrated Printer


 NEW

### Key Features

- Full body composition analysis provided in 15 seconds using precision Tanita Dual Frequency BIA Technology
- Integrated printer provides instant read out of results with topline analysis.
- Results automatically stored on the SD Card, sent to a PC or printed.
- GMon Pro Software compatible, allowing trend analysis, health risk assessments and full data management.
- Max weight capacity 270kg with 100g accuracy
- Accuracy Grade: NAWI III, MDD II-a

### Measurements

- Body fat %
  - Fat mass kg,
  - Fat free mass kg
  - Muscle mass kg
  - Total Body Water %
  - Body mass index
  - Bone mass kg
  - Physique rating
  - Visceral fat rating
  - Basal Metabolic Rate kcal
  - Basal Metabolic Rate indicator
  - Metabolic Age
- Print Out Analysis:**
    - Body Fat Analysis
    - Muscle Mass Indicator
    - BMR Indicator
    - Physique Rating
    - Target: BF and Weight

