

Chair scale KERN MCD

Mobile chair scale in ergonomically-optimised design for safe, convenient weighing



**NEW**



**Your KERN specialist dealer**

---



#### Features

- This chair scale is the ideal measuring instrument for elderly care and for paramedics, and can handle weighing overweight persons weighing more than 300 kg
- **1** Mobile version with 4 rollers and particularly convenient locking brakes on both rear wheels
- With its four wheels, this chair scale provides maximum mobility in bringing the scale to the patient. This ensures more efficient use of time for the staff and greater safety for the patients, who can be weighed in their familiar environment
- The rollers, with their large diameter, make it easier to get over door thresholds, edges and across the gaps into elevators
- For fragile patients, the comfortable, ergonomically optimised seat offers secure support during weighing
- Because of the clearly contrasting black colour of the seating surface and footrests ideally suited for use with dementia patients

- **2** Two foldable armrests and footrests make transfers into the chair easier. Ideal for overweight patients or barrier-free use, e.g. for transfers from a bed to the scale
- Ergonomically positioned carrying handles
- Hold function:  
When patients do not sit completely still, a stable weight is calculated using an average weight and this is then "frozen". This means that you have sufficient time to attend to the patient first and then take the weight reading in peace
- BMI function to determine underweight/normal weight/surplus weight
- Practical battery operation using standard batteries ensures a high level of flexibility and freedom from mains adapters, sockets, chargers, etc.
- High resolution readout: Read-out [d] for 5 sec increased by one decimal place at the touch of a key

#### Technical data

- **3** LCD display, digit height 25 mm
- Dimensions of display device W×D×H 200×130×60 mm
- Overall dimensions W×D×H 625×990×985 mm
- Batteries included, 6×1,5V AA, operating time up to 20 h
- Net weight approx. 20 kg

#### Accessories

- Protective working cover over the display device, standard, can be re-ordered, scope of delivery: 5 items, KERN MBC-A06S05
- Practical mains adapter pouch on the back of the chair as well as mains adapter (external), KERN MCC-A01
- Cleaning cloths, alcohol-free cloths for wipe disinfectant. Fulfill the legal requirements for occupational safety in accordance with TRGS 525/540. Packaging contents 80 pcs., size 20×22 cm per cloth, KERN MYC-01

STANDARD



Model	Weighing range	Readout	
	[Max] kg	[d] kg	
KERN MCD 300K-1	300	0,1	

#### Pictograms



##### Data interface RS-232:

To connect the balance to a printer, PC or network.



##### Battery operation:

Ready for battery operation. The battery type is specified for each device.



##### Weighing principle: Strain gauges

Electrical resistor on an elastic deforming body.



##### Pallet shipment:

The time required for internal shipping preparations is shown in days in the pictogram.