DC 430 MA P

Dual Frequency Body Composition Monitor with Integrated Printer

Key Features

- Full body composition analysis provided in 15 seconds using precision Tanita Dual Frequency BIA Technology
- Integrated printer provides instant read out of results with topline analysis.
- Results automatically stored on the SD Card, sent to a PC or printed.
- GMon Pro Software compatible, allowing trend analysis, health risk assessments and full data management.
- Max weight capacity 270kg with 100g accuracy
- Accuracy Grade: NAWI III, MDD II-a

Measurements

- Body fat %
- Fat mass kg,
- Fat free mass kg
- Muscle mass ka
- Total Body Water %
- Body mass index
- Bone mass kg
- Physique rating
- Visceral fat rating
- Basal Metabolic Rate kcal
- Basal Metabolic Rate indicator
- Metabolic Age

O Print Out Analysis:

- Body Fat Analysis
- Muscle Mass Indicator
- BMR Indicator
- Physique Rating
- Target: BF and Weight









WORLD No.1 BIA BRAND 5 YEAR GUARANTEE