

Crush a dose in seconds — without exposing your staff to medication dust



EasyCrush 2026 is a closed-system electric pill crusher for care homes and hospitals. Tablets are crushed directly inside a disposable medicine cup, so there is no airborne drug dust, no cross-contamination between residents, and no cleaning between doses on the medication round.

Closed system

8-second cycle

Cordless

Why care facilities choose EasyCrush

Protect your staff

Crushing tablets in an open mortar releases drug particles into the air. The closed EasyCrush cup keeps the powder contained — protecting the person preparing the dose and the next resident.

Save time on every round

Load the tablets, close the lid, and crushing starts by itself. Fit a lid and move on — no washing, no scraping, no residue carried to the next resident.

Consistent & compliant

A uniform fine powder every time, 8 to 14 selectable grinding cycles for harder tablets, and full documentation — Machinery, LVD, EMC, RoHS and WEEE — for your quality system.

Specifications

Grinding cycles	8 / 10 / 12 / 14 via MODE · standard 10 (~8 s)
Capacity	1-4 tablets (small cup) · 1-8 (large cup)
Cup sizes	Small (~30 ml) · large (~47 ml)
Controls	ON/OFF and MODE · automatic start
Dimensions	177 × 114 × 183.5 mm · fits a trolley drawer
Weight	2.1 kg
Noise level	< 48 dB
Classification	Not a medical device · HS 8438.80.9900
Power	Rechargeable battery · 4 h charge
Warranty	2 years device / 1 year battery
Country of origin	Denmark
GTIN	5745001384005
Conformity	2006/42/EC · 2014/35/EU · 2014/30/EU · 2011/65/EU · 2012/19/EU

Included & consumables

- **In the box:** device, charger, 2 grinding heads (small & large), 20 cups (10 small / 10 large), 20 rondels.
- **Medicine cups** — small & large, disposable, 2,100 pcs per box.
- **Cup lids** — LDPE, CE-marked, 2,000 pcs per box. Prepare doses ahead and transport them safely.
- **Rondels** — PP protective sheets, Ø10 cm, 2,100 pcs. Keep the grinding head clean between doses.
- **Grinding heads** — stainless steel AISI 304, Ø4.6 cm and Ø5.7 cm.

Shipping cartons (cm): device 20×20×20 (1 pc) · master 44×41×22 (4 pcs) · cups small 36×19.5×25 · cups large 47×24×30 · lids small 28.5×21.5×26 · lids large 28×26×31 · rondels 10×10×12.

Do not crush more than 8 tablets at a time. Always check the medication's own documentation, or ask a pharmacist, before crushing any tablet.

Request a quote or a 5-minute demo

Tell us how many wards or residents you supply and we reply with a quote — normally within one business day.

info@easycrush.eu

+45 52 17 58 26 · easycrush.eu

EasyCrush ApS · Nydamsvej 43 · 8362 Hørning · Denmark