## **KERN MBD 20K-2S05**



Ergonomically-optimised baby scale with large lying surface



Category	
Brand	KERN
Product categoriy	Medical scale
Product group	Baby scale
Product family	MBD

Measuring System	
Construction type of the scale	Single-range balance
Weighing system	Strain gauge
Weighing capacity [Max]	15 kg
Readability [d]	0,01 kg
Resolution	1.500
Linearity	± 0,05 kg
Reproducibility	0,01 kg
Tare range	15 kg
Default unit	kg
Units	kg Ib
Warm-up time	10 min
Eccentric loading at 1/3 [Max]	0,05 kg
Maximum creep (15 minutes)	20 g
Maximum creep (30 minutes)	40 g

	. 0
Default unit	kg
Units	kg
	lb
Warm-up time	10 min
Eccentric loading at 1/3 [Max]	0,05 kg
Maximum creep (15 minutes)	20 g
Maximum creep (30 minutes)	40 g
Approval	
CE mark	/
OL Mark	<b>Y</b>
Display	
Display Display type	LCD
	LCD ✓
Display type	
Display type Display backlight	✓
Display type Display backlight Display digit height	<b>√</b> 30 mm
Display type Display backlight Display digit height Languages of the user interface	<b>√</b> 30 mm

 $(W \times D \times H)$ 

Dimensions completely mounted (W×D×H)	597×350×74 mm
Material display housing	Plastics
Material housing	plastic
Material weighing plate	plastic
Dimensions weighing surface (W×D)	560×225 mm

- J.	•
Dimensions weighing surface (W×D)	560×225 mm
(W^D)	
Functions	
Tare function	manual (single & multi)
Power Supply	
Supplied power supply	Battery
Battery	2×1.5 V AAA
Battery / accumulator type	Alkali(-Manganese)
Battery capacity	420 mAh
Battery voltage	1,5 V
Battery operating time	20 h
Environmental conditions	
Humity of environment [Max]	80 %
Ambient temperature [Min]	10 °C
Ambient temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C
Packing & Shipping	
Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	480×654×885 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	30 kg
Shipping weight	55,6 kg
Services (optional)	
Article number for DAkkS calibration	963-128\$05
Article number for certificate of conformity (verification)	969-517

conformity (verification)

## **KERN MBD 20K-2S05**



Ergonomically-optimised baby scale with large lying surface

## Pictograms STANDARD UNIT BATT DMS 1 DAY OPTION

