

EMR-validated handrail scale with ID-Display and optional height measurement

seca Scale-up Line

- + Display shows user and patient names or IDs for increased patient safety
- + Data transmission controlled by touchscreen
- + Safe weighing with a stable handrail
- + Lightweight, sophisticated design and extremely robust thanks to steel handrail that is directly attached to the cast iron pedestal
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + Convenient transport with extremely robust rubber wheels
- + Optional panda bear facilitates measuring and weighing of children
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer



seca Scale-up Line

EMR-validated handrail scale with ID-Display and optional height measurement



Firm hold provides further stability

The elegantly-curved handrail is integrated with the cast-iron base and is part of the weight measurement, while providing added safety for patients who need it. This allows patients with restricted mobility to step onto the scale easier, the risk of falling is considerably reduced, and improves user workflow. The handrail is ergonomically designed to offer the ideal grip for patients, no matter which height.

Approval class	III
Capacity	300 kg
Graduation	50 g < 150 kg > 100 g
Dimensions (WxHxD)	801 × 1,280 × 653 mm / 22.6 × 50.4 × 25.7 in
Platform (WxHxD)	600 × 70 × 467 mm / 23.6 × 2.8 × 18.4 in
Net weight	25 kg / 55.1 lbs
Power supply	Power adapter
Functions	Display of user and patient name or ID, Display date of birth patient, TARE, Pre-TARE, HOLD, Auto-HOLD, CLEAR, Auto-CLEAR, RESET, Automatic weighing range switch-over, Adjustable damping, BMI, Mother+child function, Standby mode
Interfaces	Ethernet, USB for barcode scanner, Wi-Fi
Accessories	+ Ultrasonic height measurement seca 257 + Panda seca 459 + Barcode scanner mount seca 463
Medical Device Regulation (MDR) certified	



Spacious and hygienic platform

The platform of the seca Scale-up Line consists of absolutely break-proof seca Bearclaw Glass®. The glass not only looks good, it also ages significantly slower than e.g. plastic and is particularly hygienic. It can be easily and quickly cleaned and has germ-reduction rates similar to stainless steel, with an anti-slip coating for safe standing. Transport is very easy due to the smooth-running, robust rubber wheels.



ID-Display ensures easy and user-friendly data transmission to EMR systems

seca Scale-up Line is developed especially for perfect integration of measured data into the EMR. Weight and data such as time of measurement are stored centrally and accessible for every user. The ID-Display shows name or patient ID with every measurement, which are gathered by an additional barcode scanner. This makes the integration process transparent, safe and even more user-friendly.

Order code **ON SHM ANT NN**

Flat scale seca 655
+ Handrail seca 455